

WESAD Exploration

March 6, 2025

```
[186]: import pandas as pd
import matplotlib.pyplot as plt
import numpy as np
import flirt
import neurokit2 as nk

%matplotlib widget
```

Load wesad data specimen.

```
[187]: wesad_data = '/home/art/Downloads/WESAD/S7/S7.pkl'
dataset = pd.read_pickle(wesad_data)
BVP = dataset['signal']['wrist']['BVP'][:,0]
EDA = dataset['signal']['wrist']['EDA'][:,0]
TEMP = dataset['signal']['wrist']['TEMP'][:,0]
RESP = dataset['signal']['chest']['Resp'][:,0]
ECG = dataset['signal']['chest']['ECG'][:,0]
label = dataset['label']
```

Upsample wrist data.

```
[188]: ECG_down = nk.signal_resample(ECG, sampling_rate=700,
    ↳desired_sampling_rate=256, method='interpolation')
RESP_down = nk.signal_resample(RESP, sampling_rate=700,
    ↳desired_length=len(ECG_down), method='interpolation')
BVP_up = nk.signal_resample(BVP, desired_length=len(ECG_down), method='FFT')
EDA_up = nk.signal_resample(EDA, desired_length=len(ECG_down), method='poly')
TEMP_up = nk.signal_resample(TEMP, desired_length=len(ECG_down), method='poly')
label_resampled = nk.signal_resample(label, desired_length=len(ECG_down),
    ↳method='interpolation')

data_df = pd.DataFrame(
    {
        'BVP': BVP_up,
        'EDA': EDA_up,
        'TEMP': TEMP_up,
        'ECG': ECG_down,
        'RESP': RESP_down,
```

```

        'label' : label_resampled,
    }
)

data_df.head()

```

```

[188]:
      BVP      EDA      TEMP      ECG      RESP  label
0  9.570000  5.973745  33.272443 -0.019226  1.832581      0
1  13.690549  6.038064  33.622649 -0.012980  1.837982      0
2  13.802470  6.099970  33.959122 -0.031913  1.795855      0
3  11.086804  6.159436  34.281717 -0.110225  1.797355      0
4   7.410000  6.216439  34.590313 -0.095998  1.817480      0

```

Get Features out of WESAD Data.

```

[189]:
ecg_signals, ecg_info = nk.ecg_process(data_df['ECG'], sampling_rate=256)
nk.ecg_plot(ecg_signals, ecg_info)
eda_signals, eda_info = nk.eda_process(data_df['EDA'], sampling_rate=256,
    ↪report='Text')
nk.eda_plot(eda_signals, eda_info)
ppg_signals, ppg_info = nk.ppg_process(data_df['BVP'], sampling_rate=256,
    ↪report='Text')
nk.ppg_plot(ppg_signals, ppg_info)
resp_signals, resp_info = nk.rsp_process(data_df['RESP'], sampling_rate=256,
    ↪method='khodadad2018', method_rvt='harrison2021', report='Text')
nk.rsp_plot(resp_signals, resp_info)

```

The raw signal, sampled at 256 Hz, was cleaned using the default method of the neurokit2 package.

The signal was decomposed into phasic and tonic components using the method described in neurokit.

The cleaned signal was used to detect peaks using the default method of the `neurokit2` package.

References

- <https://doi.org/10.21105/joss.01667>

The raw signal, sampled at 256 Hz, was preprocessed using a bandpass filter ([0.5 - 8 Hz], Butterworth 3rd order; following Elgendi et al., 2013).

The peak detection was carried out using the method described in Elgendi et al. (2013).

The quality assessment was carried out using template-matching, approximately as described in Orphanidou et al. (2015).

	PPG_Rate_Mean	PPG_Rate_SD
0	66.285141	24.456355

References

- Elgendi M, Norton I, Brearley M, Abbott D, Schuurmans D (2013)
Systolic Peak Detection in Acceleration Photoplethysmograms
Measured from Emergency Responders in Tropical Conditions
PLoS ONE 8(10): e76585. doi:10.1371/journal.pone.0076585.
- Makowski, D., Pham, T., Lau, Z. J., Brammer, J. C., Lespinasse, F., Pham, H.,
Schölzel, C., & Chen, S. A. (2021). NeuroKit2: A Python toolbox for
neurophysiological signal processing.
Behavior Research Methods, 53(4), 1689-1696.
<https://doi.org/10.3758/s13428-020-01516-y>
- Orphanidou C, Bonnici T, Charlton P, Clifton D, Vallance D, Tarassenko L
(2015)
Signal-quality indices for the electrocardiogram and
photoplethysmogram: Derivation
and applications to wireless monitoring
IEEE Journal of Biomedical and Health Informatics 19(3): 832-838.
doi:10.1109/JBHI.2014.2338351.
The raw signal, sampled at 256 Hz, was preprocessed using a second order 0.05-3
Hz bandpass Butterworth filter.

The peak detection was carried out using the method described in Khoadadad et al. (2018).

The respiratory volume per time computation was carried out using the method described in Harrison et al. (2021).


	RSP_Rate_Mean	RSP_Rate_SD
0	15.502424	6.14461

References

- Harrison, S. J., Bianchi, S., Heinzle, J., Stephan, K. E., Iglesias, S., &
Kasper, L. (2021).
A Hilbert-based method for processing respiratory timeseries.
Neuroimage, 230, 117787.
- Khodadad, D., Nordebo, S., Müller, B., Waldmann, A., Yerworth, R., Becher, T.,
... & Bayford, R. (2018).
Optimized breath detection algorithm in electrical impedance
tomography.
Physiological measurement, 39(9), 094001.
- Makowski, D., Pham, T., Lau, Z. J., Brammer, J. C., Lespinasse, F., Pham, H.,
Schölzel, C., & Chen, S. A. (2021). NeuroKit2: A Python toolbox for
neurophysiological signal processing.
Behavior Research Methods, 53(4), 1689-1696.


Electrocardiogram (ECG)

ECG signal and peaks individual Heart Beats (average heart rate:




Electrodermal Activity (EDA)


Raw and Cleaned Signal



Skin Conductance Response (SCR)





Skin Conductance Level (SCL)



Photoplethysmogram (PPG)

PPG signal and individual Heart Beats (average heart rate:





```
[190]: quick_feat = flirt.with_.empatica('/home/art/Downloads/WESAD/S7/S7_E4_Data.zip',
                                         window_length= 60,
                                         window_step_size = 1,
                                         hrv_features = True,
                                         eda_features = True,
                                         acc_features = True)
```

```
HRV features: 100%|          | 6403/6403 [00:05<00:00, 1116.35it/s]
```

```
EDA features:  0%|          | 0/6494 [00:00<?, ?it/s]
```

```
/tmp/ipykernel_335021/1998099736.py:1: UserWarning: phasic_entropy contains more
than 5% (actual: 63.3%) nan, inf, or -inf values. We recommend to delete this
feature column.
```

```
quick_feat =
flirt.with_.empatica('/home/art/Downloads/WESAD/S7/S7_E4_Data.zip',
```

```
ACC features:  0%|          | 0/6495 [00:00<?, ?it/s]
```

```
/home/art/miniconda3/envs/biosig/lib/python3.12/site-
packages/flirt/stats/common.py:35: RuntimeWarning: Precision loss occurred in
moment calculation due to catastrophic cancellation. This occurs when the data
are nearly identical. Results may be unreliable.
```

```

    results[key] = value(data)
/home/art/miniconda3/envs/biosig/lib/python3.12/site-
packages/flirt/stats/common.py:35: RuntimeWarning: Precision loss occurred in
moment calculation due to catastrophic cancellation. This occurs when the data
are nearly identical. Results may be unreliable.
    results[key] = value(data)

```

```

[192]: from skchange.change_detectors import MovingWindow
       from skchange.change_scores import CUSUM

       detector = MovingWindow(
           change_score=CUSUM(), # measures the difference in mean between two
           ↪ intervals.
           bandwidth=20,
           threshold_scale=1.0,
       )
       display(detector)

```

```
MovingWindow(bandwidth=20, change_score=CUSUM(), threshold_scale=1.0)
```

```


[193]: hrv_hf = quick_feat['hrv_hf'].dropna().reset_index(drop=True)
       hrv_hf_detections = detector.fit_predict(hrv_hf)
       hrv_changepoints = hrv_hf_detections.get("ilocs", [])

       fig, ax = plt.subplots(figsize=(10, 5))
       ax.plot(hrv_hf.index, hrv_hf.values, label="HRV_HF Signal", color="blue")

       for i, hrv_changepoint in enumerate(hrv_changepoints):
           ax.axvline(x=hrv_changepoint, linestyle="dashed", color="red",
           ↪ label="Changepoint" if i == 0 else "")

       ax.set_title("HRV_HF Signal with Detected Changepoints")
       ax.set_xlabel("Sample Index")
       ax.set_ylabel("HRV_HF Value")
       ax.legend()
       plt.show()

```


```
[194]: detector = MovingWindow(change_score=CUSUM(), bandwidth=180, threshold_scale=1.0)
display(detector)
eda_phasic_mean = quick_feat['eda_phasic_mean'].dropna().reset_index(drop=True)
eda_phasic_detect = detector.fit_predict(eda_phasic_mean)
eda_phasic_changepoints = eda_phasic_detect.get("ilocs", [])

fig, ax = plt.subplots(figsize=(10, 5))
ax.plot(eda_phasic_mean.index, eda_phasic_mean.values, label="EDA Phasic Mean",
        color="blue")

for i, eda_phasic_changepoints in enumerate(eda_phasic_changepoints):
    ax.axvline(x=eda_phasic_changepoints, linestyle="dashed", color="red",
               label="Changepoint" if i == 0 else "")

ax.set_title("EDA Phasic with Detected Changepoints")
ax.set_xlabel("Sample Index")
ax.set_ylabel("HRV entropy Value")
ax.legend()
plt.show()
```

```
MovingWindow(bandwidth=180, change_score=CUSUM(), threshold_scale=1.0)
```




```
[195]: detector = MovingWindow(change_score=CUSUM(), bandwidth=180, threshold_scale=1)
display(detector)
hrv_energy = quick_feat['hrv_energy'].dropna().reset_index(drop=True)
hrv_energy_detect = detector.fit_predict(hrv_energy)
hrv_energy_changepoints = hrv_energy_detect.get("ilocs", [])

fig, ax = plt.subplots(figsize=(10, 5))
ax.plot(hrv_energy.index, hrv_energy.values, label="HRV Energy", color="blue")

for i, hrv_energy_changepoints in enumerate(hrv_energy_changepoints):
    ax.axvline(x=hrv_energy_changepoints, linestyle="dashed", color="red",
    label="Changepoint" if i == 0 else "")

ax.set_title("HRV Energy with Detected Changepoints")
ax.set_xlabel("Sample Index")
ax.set_ylabel("HRV Energy Value")
ax.legend()
plt.show()
```

MovingWindow(bandwidth=180, change_score=CUSUM(), threshold_scale=1)




```
[196]: detector = MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=1.
    ↪0)
display(detector)
ECG_IBIs = quick_feat['num_ibis'].dropna().reset_index(drop=True)
ecg_ibi_detect = detector.fit_predict(ECG_IBIs)
ecg_ibi_changepoints = ecg_ibi_detect.get("ilocs", [])

fig, ax = plt.subplots(figsize=(10, 5))
ax.plot(ECG_IBIs.index, ECG_IBIs.values, label="HRV IBIs", color="blue")

for i, ecg_ibi_changepoints in enumerate(ecg_ibi_changepoints):
    ax.axvline(x=ecg_ibi_changepoints, linestyle="dashed", color="red",
    ↪label="Changepoint" if i == 0 else "")

ax.set_title("ECG IBIs with Detected Changepoints")
ax.set_xlabel("Sample Index")
ax.set_ylabel("ECG IBI Value")
ax.legend()
plt.show()
```

```
MovingWindow(bandwidth=1000, change_score=CUSUM(), threshold_scale=1.0)
```




```
[197]: detector = MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=1.
      ↪0)
display(detector)
hr_max = quick_feat['hrv_max_hr'].dropna().reset_index(drop=True)
hr_max_detect = detector.fit_predict(hr_max)
hr_max_changepoints = hr_max_detect.get("ilocs", [])

fig, ax = plt.subplots(figsize=(10, 5))
ax.plot(hr_max.index, hr_max.values, label="HR Max", color="blue")

for i, hr_max_changepoints in enumerate(hr_max_changepoints):
    ax.axvline(x=hr_max_changepoints, linestyle="dashed", color="red",
    ↪label="Changepoint" if i == 0 else "")

ax.set_title("HR Max with Detected Changepoints")
ax.set_xlabel("Sample Index")
ax.set_ylabel("HR Max Value")
ax.legend()
plt.show()
```

```
MovingWindow(bandwidth=1000, change_score=CUSUM(), threshold_scale=1.0)
```



```
[ ]: fig, axes = plt.subplots(nrows=7, ncols=1, figsize=(12, 22), sharex=True)
fig.suptitle("Subject 7")
signals = [
    ("hrv_max_hr", "HR Max"),
    ("num_ibis", "ECG IBIs"),
    ("hrv_energy", "HRV Energy"),
    ("eda_phasic_mean", "EDA Phasic Mean"),
    ("hrv_hf", "HRV_HF Signal"),
    ("hrv_rmssd", "HRV RMSSD"),
    ("hrv_lf_hf_ratio", "HRV LF/HF Ratio")
]
# Hyperparameters....
detectors = [
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=1.0),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=1.0),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=1.0),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=1.0),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=1.0),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=1.0),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=1.0)
]

for ax, (signal_name, title), detector in zip(axes, signals, detectors):
    signal_data = quick_feat[signal_name].dropna().reset_index(drop=True)
    detection_results = detector.fit_predict(signal_data)
    changepoints = detection_results.get("ilocs", [])

    ax.plot(signal_data.index, signal_data.values, label=title, color="blue")
```

```

    for i, changepoint in enumerate(changepoints):
        ax.axvline(x=changepoint, linestyle="dashed", color="red",
        ↪label="Changepoint" if i == 0 else "")

    ax.set_title(f"{title} with Detected Changepoints")
    ax.set_ylabel(f"{title} Value")
    ax.legend()

axes[-1].set_xlabel("Sample Index")

plt.show()

```

Subject 4

```

[199]: quick_feat4 = flirt.with_.empatica('/home/art/Downloads/WESAD/S4/S4_E4_Data.
        ↪zip',

                                window_length= 60,
                                window_step_size = 1,
                                hrv_features = True,
                                eda_features = True,
                                acc_features = True,
                                debug = True
                                )

```

Reading files

Calculating HRV features

HRV features: 100% | 7771/7771 [00:00<00:00, 12245.69it/s]

Calculating EDA features

EDA features: 0% | 0/8000 [00:00<?, ?it/s]

Calculating ACC features

/tmp/ipykernel_335021/419229283.py:1: UserWarning: phasic_entropy contains more than 5% (actual: 73.44%) nan, inf, or -inf values. We recommend to delete this feature column.

```

    quick_feat4 =
    flirt.with_.empatica('/home/art/Downloads/WESAD/S4/S4_E4_Data.zip',

```

ACC features: 0% | 0/8000 [00:00<?, ?it/s]

/home/art/miniconda3/envs/biosig/lib/python3.12/site-packages/flirt/stats/common.py:35: RuntimeWarning: Precision loss occurred in moment calculation due to catastrophic cancellation. This occurs when the data are nearly identical. Results may be unreliable.

```

    results[key] = value(data)

```

```

[200]: fig, axes = plt.subplots(nrows=7, ncols=1, figsize=(12, 20), sharex=True)
fig.suptitle('Subject 4')
signals = [
    ("hrv_max_hr", "HR Max"),
    ("num_ibis", "ECG IBIs"),
    ("hrv_energy", "HRV Energy"),
    ("eda_phasic_mean", "EDA Phasic Mean"),
    ("hrv_hf", "HRV_HF Signal"),
    ("hrv_rmssd", "HRV RMSSD"),
    ("hrv_lf_hf_ratio", "HRV LF/HF Ratio")
]
# Hyperparameters....
detectors = [
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=.5),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=.5),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=.5),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=.5),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=.5),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=.5),
    MovingWindow(change_score=CUSUM(), bandwidth=1000, threshold_scale=.5)
]

for ax, (signal_name, title), detector in zip(axes, signals, detectors):
    signal_data = quick_feat4[signal_name].dropna().reset_index(drop=True)
    detection_results = detector.fit_predict(signal_data)
    changepoints = detection_results.get("ilocs", [])

    ax.plot(signal_data.index, signal_data.values, label=title, color="blue")


    for i, changepoint in enumerate(changepoints):
        ax.axvline(x=changepoint, linestyle="dashed", color="red",
        label="Changepoint" if i == 0 else "")

    ax.set_title(f"{title} with Detected Changepoints")
    ax.set_ylabel(f"{title} Value")
    ax.legend()

axes[-1].set_xlabel("Sample Index")
plt.show()

```

Subject 4



Subject 6

```
[201]: quick_feat6 = flirt.with_.empatica('/home/art/Downloads/WESAD/S6/S6_E4_Data.
      ↪zip',
      window_length= 60,
      window_step_size = 1,
      hrv_features = True,
      eda_features = True,
      acc_features = True,
      debug = True
      )
```

Reading files

Calculating HRV features

HRV features: 100%| | 8085/8085 [00:00<00:00, 40700.97it/s]

Calculating EDA features

EDA features: 0%| | 0/8324 [00:00<?, ?it/s]

Calculating ACC features

/tmp/ipykernel_335021/1399992444.py:1: UserWarning: phasic_entropy contains more than 5% (actual: 63.27%) nan, inf, or -inf values. We recommend to delete this feature column.

```
    quick_feat6 =
    flirt.with_.empatica('/home/art/Downloads/WESAD/S6/S6_E4_Data.zip',
```

ACC features: 0%| | 0/8324 [00:00<?, ?it/s]

```
[202]: fig, axes = plt.subplots(nrows=7, ncols=1, figsize=(12, 20), sharex=True)
fig.suptitle('Subject 6')
signals = [
    ("hrv_max_hr", "HR Max"),
    ("num_ibis", "ECG IBIs"),
    ("hrv_energy", "HRV Energy"),
    ("eda_phasic_mean", "EDA Phasic Mean"),
    ("hrv_hf", "HRV_HF Signal"),
    ("hrv_rmssd", "HRV RMMSD"),
    ("hrv_lf_hf_ratio", "HRV LF/HF Ratio")
]
# Hyperparamaters....
detectors = [
    MovingWindow(change_score=CUSUM(), bandwidth=1800, threshold_scale=.5),
    MovingWindow(change_score=CUSUM(), bandwidth=1800, threshold_scale=.5),
    MovingWindow(change_score=CUSUM(), bandwidth=1800, threshold_scale=.5),
    MovingWindow(change_score=CUSUM(), bandwidth=1800, threshold_scale=.5),
```

```

MovingWindow(change_score=CUSUM(), bandwidth=1800, threshold_scale=.5),
MovingWindow(change_score=CUSUM(), bandwidth=1800, threshold_scale=.5),
MovingWindow(change_score=CUSUM(), bandwidth=1800, threshold_scale=.5)
]

for ax, (signal_name, title), detector in zip(axes, signals, detectors):
    signal_data = quick_feat6[signal_name].dropna().reset_index(drop=True)
    detection_results = detector.fit_predict(signal_data)
    changepoints = detection_results.get("ilocs", [])

    ax.plot(signal_data.index, signal_data.values, label=title, color="blue")

    for i, changepoint in enumerate(changepoints):
        ax.axvline(x=changepoint, linestyle="dashed", color="red",
↪label="Changepoint" if i == 0 else "")

    ax.set_title(f"{title} with Detected Changepoints")
    ax.set_ylabel(f"{title} Value")
    ax.legend()

axes[-1].set_xlabel("Sample Index")
plt.show()

```

Subject 6

